RS VIVO Model Natural Pink

INSTRUCTIONS FOR USE

rapidshape

PACKAGING UNITS

RS VIVO Model Natural Pink, 1000g [RS006191]

DESCRIPTION

RS VIVO Model Natural Pink is a flowable, light-curing, acrylate-based composite material for the additive manufacturing of dental models. **RS VIVO Model Natural Pink** is intended to be cured by exposure to LED light. The resin has been optimized for rapidshape® 3D printing systems at wavelengths of 385/405 nm for layer thicknesses of 50/100 µm.

COMPATIBLE EQUIPMENT

RS VIVO Model Natural Pink is compatible with the following equipment:

3D printer: rapidshape® DII, D+, PRO, ONE Series 3D printing systems with wavelengths of 385/405 nm.

Washing unit: rapidshape® RS wash, WASH with appropriate parameter set.

Post curing unit: rapidshape® RS cure, RS cure XL, CURE with appropriate parameter set; Otofash G171 Cure Box/NK Optik.

The user must follow all applicable labeling including the instructions for use, installation guides, and other related labeling for all listed compatible devices used in conjunction with **RS VIVO Model Natural Pink**. Strict adherence to all applicable instructions is essential to ensure safe printing results.

Ensure that equipment is installed, maintained, and calibrated according to the manufacturer's instructions for use. Use only validated and approved equipment as listed above.

COMPOSITION

Composed of acrylates and initiators. Always refer to the latest safety data sheet!

INDICATION

RS VIVO Model Natural Pink is a light-curing dental resin for the additive manufacturing of dental models.

ADVERSE REACTIONS

User: Product may cause allergic reactions. Observe warnings.

APPLICATION AND CONSTRUCTION PROCESS

PREPARE DATA

In the slicing software, place and orient the parts on the build platform. Using predefined workflow parameters, parts are automatically placed on the platform in the best possible way without any additional steps necessary.

SUPPORT CONSTRUCTION

For best results, it is recommended that all parts be supported with a bar contour cross, and a hexagonal grid as a baseplate. Using predefined workflow parameters, parts are automatically supported in the best possible way without any additional steps necessary. Do not place supports on the preparation surfaces as they will interfere with the fit.

ESTABLISH VARIABLE PROCESS PARAMETERS

Select the appropriate process printing parameters according to the system specifications. Ensure that the correct material-specific parameter set is selected. Ensure that only product-specific, pre-validated settings are used.

JOB PREVIEW

Examine each layer for problem features, especially overhangs that need support. Find and correct areas that are not effectively connected to the build platform and/or the intended supports.

MATERIAL USE

General notes:

- Make sure the resin is at an ambient temperature (20-25°C/68-77°F) before starting printing.
- Do not use the resin after its expiration date.
- Clean the build platform after each job.
- In the event of a print failure, filter the resin through a fine sieve with a pore size of at least 190 µm.
- It is recommended to filter the resin frequently.

INITIAL FILLING

- Be sure to use a clean material reservoir.
- Stir the resin in the bottle thoroughly with a spatula before filling it into the reservoir. Avoid bubbles and splashes.
- Allow the resin to stand for 10 minutes after filling before starting the print job.

REUSE OF MATERIAL

- Ensure that the resin in the material reservoir as well as the reservoir bottom is free of residues from previous print jobs. If necessary, use paper cards/a rubber spatula to remove cured particles or use the «Automatic Reservoir Cleaning»/«Clean Reservoir» function of the printer (refer to the printer manufacturer's instructions for use).
- Mix the resin thoroughly with a paper card/a rubber spatula in the reservoir. Avoid bubbles and splashes.
- Allow the resin to stand for 10 minutes after mixing before starting the print job.

REFILL OF MATERIAL

- Ensure that the resin in the material reservoir as well as the reservoir bottom is free of residues from previous print jobs. If necessary, use paper cards/a rubber spatula to remove cured particles or use the «Automatic Reservoir Cleaning»/«Clean Reservoir» function of the printer (refer to the printer manufacturer's instructions for use).
- Stir the resin in the bottle thoroughly with a spatula and the remaining resin in the material reservoir with a paper card/a rubber spatula.
- Add the prepared resin to the reservoir. Avoid bubbles and splashes.
- Allow the resin to stand for 10 minutes after filling before starting the print job.

START PRNT JOB

CLEANING AND DRYING

Option 1:

Clean the built part under flowing conditions using the washing unit rapidshape® RS wash or WASH with the appropriate parameter set.

Option 2:

Clean the built part under flowing conditions with isopropanol for a maximum of 5 minutes. If necessary, carefully clean the part with a brush and a spray bottle filled with isopropanol.

Blow the part with compressed air and allow it to dry for 30 minutes at room temperature. Alternatively, the part can be stored in a drying cabinet at a temperature of approximately 37°C/99°F.

POST CURING

For post curing, the built part must be dry and cooled to room temperature.

Option 1:

Cure the built part using the curing unit rapidshape® RS cure, RS cure XL or CURE with the appropriate parameter set.

Option 2:

Cure the built part using the Otofash G171 Cure Box with 2000 flashes each for top and bottom (Σ4000). (Recommendation: Curing under exclusion of oxygen prevents the formation of the inhibition layer.)

Ensure that all equipment used is maintained and calibrated in accordance with the manufacturer's instructions for use.

FINAL DRYING

Blow the part with compressed air and allow it to dry for 30 minutes. Alternatively, the part can be stored in a drying cabinet at a temperature of approximately 37°C/99°F.

REMOVE SUPPORT STRUCTURES

Remove supports and grind remnants to achieve a smooth surface.

FURTHER PROCESSING INSTRUCTIONS

The **RS VIVO Model Natural Pink** restorations are subsequently processed and finished according to the normal dental technical processes.

WARNINGS

Applies to the uncured material: May cause allergic skin reactions. Avoid breathing vapors. Ensure adequate ventilation. Avoid contact with uncured material. Wear protective gloves/protective clothing/eye protection. If on skin: Wash with plenty of water and soap. Remove contaminated clothing and wash before reuse. If skin irritation or rash occurs: Get medical advice/attention. After contact with the eyes, rinse open eyes for several minutes under running water and seek medical advice if irritation persists. The unpolymerized material is harmful to the environment. Avoid release to the environment. Collect any spillage. Always observe the relevant country-specific health and safety regulations and legal requirements when handling hazardous substances. Always refer to the latest safety data sheet!

STORAGE INSTRUCTIONS / SHELF LIFE

Protect from exposure to light. Keep the container tightly sealed. May spontaneously polymerize. Store between 4°C and 28°C (39°F-82°F). Do not use after the expiration date. See imprint/label for expiration date.

DISPOSAL INSTRUCTIONS

Dispose of contents/container in accordance with local regulations.

IMPORTANT PRODUCT INFORMATION

The product has been developed for use in dentistry and must be used in accordance with the instructions for use.

-  Batch code
-  Catalogue number
-  Use-by date
-  Consult instructions for use
-  Temperature limit
-  Keep away from sunlight
-  Stir before use
-  Manufacturer
-  Exclamation mark
-  Environmental hazard

 **Rapid Shape GmbH**
Roemerstr. 21
71296 Heimsheim/Germany
+49 7033 309 878-0
www.rapidshape.de

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