**rs  
VIVO**

# **RS VIVO Denture Impact**

Pink  
Pink-Transparent

**rapidshape**

|          |                           |    |
|----------|---------------------------|----|
| bg       | Инструкция за употреба    | 4  |
| cs       | Návod k použití           | 7  |
| da       | Brugsanvisning            | 10 |
| de       | Gebrauchsanweisung        | 13 |
| el       | Οδηγίες χρήσης            | 16 |
| en       | Instructions for use      | 19 |
| en [USA] | Instructions for use      | 22 |
| es       | Modo de empleo            | 28 |
| et       | Kasutusjuhend             | 31 |
| fi       | Käyttöohje                | 34 |
| fr       | Mode d'emploi             | 37 |
| hr       | Upute za uporabu          | 40 |
| hu       | Használati útmutató       | 43 |
| it       | Istruzioni per l'uso      | 46 |
| lt       | Naudojimo instrukcija     | 49 |
| lv       | Lietošanas instrukcija    | 52 |
| nl       | Gebruiksaanwijzing        | 55 |
| no       | Bruksanvisning            | 58 |
| pl       | Instrukcja stosowania     | 61 |
| pt       | Manual de instruções      | 64 |
| ro       | Instrucțiuni de utilizare | 67 |
| ru       | Инструкция по применению  | 70 |
| sk       | Návod na použitie         | 73 |
| sl       | Navodila za uporabo       | 76 |
| sv       | Bruksanvisning            | 79 |
| tr       | Kullanım kılavuzu         | 82 |

## ORDERING INFORMATION

RS VIVO Denture Impact, 1000 g

Part no.

---

|      |                                        |
|------|----------------------------------------|
| Pink | 04720 [Rapid Shape part no.: RS006200] |
|------|----------------------------------------|

---

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Denture bases

**Patient target group**

Persons being treated in the context of a dental procedure.

**Intended users**

Dentist, dental technician

**Suitable for the following DLP printers/cleaning/post-exposure**

see "Annex 1" (enclosed separately)

**Processing**

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Polish surface mechanically.
- Processing temperature  $23\text{ °C} \pm 2\text{ °C}$ .

**Safety information**

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- MMA free.
- Refer to the relevant safety data sheet for hazard and safety information.

## Cleaning

### By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

### By the patient

For daily cleaning of the prosthesis, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards. Denture cleaning tablets can be used following the manufacturer's instructions.

### Note

- Do not use toothpaste; the abrasive particles are coarse and scratch the denture.
- Do not use any acids or solvents.

## Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Read and understand the safety data sheet!
- For users and/or patients

Any serious incidents occurring in relation to this product should be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent authority of the Member State in which the user and/or patient is established.

## Storage

- RS VIVO Denture Impact is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

## Contraindication

Contains (meth)acrylics and phosphine oxide. Some ingredients of RS VIVO Denture Impact may cause allergic reactions in predisposed persons. In such cases refrain from using the product. RS VIVO Denture Impact must only be inserted intraorally in completely polymerised state.

## Adverse effects

Product may cause allergic reactions.

## Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

## Storage

15 °C  
59 °F

28 °C  
82 °F

Processing  
at 23 °C ± 2 °C

## Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer.

### ① Construction process

Generation of a Print Job complying with machine and material parameters.

### ② Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

### ③ Cleaning

see “Annex 1, Cleaning equipment”

### ④ Post-exposure

see “Annex 1, Light curing equipment”

### ⑤ Surface processing

Polish surface mechanically.

**Caution statement**

Caution: Federal Law restricts this device to sale by or on the order of a dentist.

**Indications for use**

The RS VIVO Denture Impact is a light-cured resin indicated for the fabrication of all types of denture bases, for instance full and partial removable dentures. The RS VIVO Denture Impact is intended for continuous use in the oral environment, exclusively for professional dental work. RS VIVO Denture Impact can be used in combination with a stereolithographic 3D printer using a 385 nm light source. A 3D-printer is not part of the device.

**Patient target group**

Persons being treated in the context of a dental procedure.

**Intended users**

Dentist, dental technician

**Special manufacturing requirements**

The printing machine and post-curing unit should be set-up, validated and maintained according to their labelling and instructions for use. RS VIVO Denture Impact is only intended to be used with the identified compatible equipment in Annex 1.

**Environmental conditions**

Processing temperature  $23^{\circ}\text{C} \pm 2^{\circ}\text{C} / 73^{\circ}\text{F} \pm 4^{\circ}\text{F}$ .

**Cleaning supplies**

Rinse bathtubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, scraper

**Processing**

- The use of RS VIVO Denture Impact is only approved when applied with the compatible devices mentioned in this instruction for use as outlined in Annex 1.
- Generate the object (STL-file) using a commercial CAD software, which is intended for dental applications.
- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the exposure device is in an orderly condition and that the moulds are completely cured (observe manufacturing process).
- After storage, the material in the bottle should be shaken and homogenized with a bottle roller before use.
- Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.

**Technical Data**

ISO 20795-1-Type 4

Does not contain MMA monomer.

Does not contain phthalates.

**Processing**

at  $23^{\circ}\text{C} \pm 2^{\circ}\text{C} / 73^{\circ}\text{F} \pm 4^{\circ}\text{F}$

**Storage**

$15^{\circ}\text{C}$   
59 °F

$28^{\circ}\text{C}$   
82 °F

### Safety information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- If resin gets into your eyes immediately rinse and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- Refer to the relevant safety data sheet for hazard and safety information.

### Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Read and understand the safety data sheet!

### MRI Safety Information

RS VIVO Denture Impact is MR Safe as it is composed of materials that are electrically nonconductive, nonmetallic and nonmagnetic.

### Storage

RS VIVO Denture Impact is to be stored dry at 15 °C - 28 °C / 59 °F - 82 °F and protected from light. Minimal influence of light can already induce polymerisation. Please cover the material vat with its lid or a glass plate to protect material from contamination.

### Contraindications

Contains (meth)acrylics and phosphine oxide.

Some ingredients of RS VIVO Denture Impact may cause allergic reactions in predisposed persons. In such cases refrain from using the product. RS VIVO Denture Impact only insert intraorally in completely polymerised state.

**Adverse effects**

Product may cause allergic reactions.

**Disposal**

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

**Material properties RS VIVO Denture Impact**

|                     |                                 |
|---------------------|---------------------------------|
| Color*              | pink, pink-transparent          |
| Viscosity*          | 400-1000 mPas                   |
| Flexural strength** | > 65 mPa                        |
| Flexural modulus**  | > 2000 mPa                      |
| Water sorption**    | < 32 $\mu\text{g}/\text{mm}^3$  |
| Water solubility**  | < 1,6 $\mu\text{g}/\text{mm}^3$ |

\*applies to liquid resin; \*\*applies to cured objects

## Manufacturing process

### ① Mixing

Before filling resin into the reservoir, the material must be mixed. Use a bottle roller or mix/shake it by hand. Inadequate mixing could cause deviations of color and print failures.

### ② Filling of printer

Fill 3D printing resin in the reservoir of the printer. It is not allowed to pour 3D printing resin from the reservoir of the printer back into the resin bottle.

### ③ Construction process

After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use. Generate a print job complying with printer and material parameters. For setting up the printer follow the instructions for use of the printer. Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer. Remaining print resin material after printing is not allowed to be poured back into the resin bottle.

### ④ Post-processing

If possible, post-processing should commence immediately following with this construction process. After raising the platform, a drip time of approx. 10 minutes is recommended. Remove the platform from the printer and remove the construction components using appropriate device (e.g., knife or spatula). Carefully remove excess resin using a light flow of compressed air.

### ⑤ Cleaning

The construction components should be cleaned in two steps with isopropyl alcohol (purity  $\geq 98\%$ ) using an ultrasonic bath. For additional suitable cleaning unit refer to Annex 1, item "Cleaning equipment".

1. Clean the construction components for 3 minutes in a reusable isopropyl alcohol (purity  $\geq 98\%$ ) using an ultrasonic bath. Then carefully remove the construction components from the support structure.
2. The precleaned construction components must be cleaned thoroughly for 3 minutes using a fresh isopropyl alcohol (purity  $\geq 98\%$ ) with an ultrasonic bath. Prior to post-exposure, check the construction components for residues. Then blow off with compressed air. After cleaning the surface of the construction components must no longer be sticky (dry) and not shiny (matt).

### ⑥ Post-exposure

Post-exposure is performed with a suitable curing unit (See Annex 1, item „Light curing equipment“).

## Caution

The use of RS VIVO Denture Impact is only approved when applied with the compatible devices mentioned in this instruction for use. Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specification.

## Finishing

Finish the construction components by using conventional dental methods and instruments.

## Symbol explanation

| Graphic                                                                         | Reference number | Title                                                                   | Description                                                                                 | Reference                                                                                                                                   |
|---------------------------------------------------------------------------------|------------------|-------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
|  | 5.1.1            | Manufacturer                                                            | Indicates the medical device manufacturer.                                                  | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.1.4            | Use-by date                                                             | Indicates the date after which the medical device is not to be used.                        | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.1.5            | Batch code                                                              | Indicates the manufacturer's batch code so that the batch or lot can be identified.         | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.1.6            | Catalogue number                                                        | Indicates the manufacturer's catalogue number so that the medical device can be identified. | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.1.9            | Distributor                                                             | Indicates the entity distributing the medical device into the locale.                       | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.3.2            | Keep away from sunlight                                                 | Indicates a medical device that needs protection from light sources.                        | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.3.7            | Temperature limit                                                       | Indicates the temperature limits to which the medical device can be safely exposed.         | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | 5.4.3            | Consult instructions for use or consult electronic instructions for use | Indicates the need for the user to consult the instructions for use.                        | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |

| Graphic                                                                          | Title              | Description                                                                                                                                                                              | Reference                                                                                                                             |
|----------------------------------------------------------------------------------|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
|   | Prescription only  | Requires a prescription in the United States of America<br>Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of, a dentist. | USA Code of Federal Regulations<br>21 CFR Part 801<br>§ 801.109 (b)(1)                                                                |
|   | Magnetic resonance | The RS VIVO Denture Impact is MR Safe.                                                                                                                                                   | FDA Guidance document: Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment dated May 20, 2021. |
|  | Warning            | The uncured material is considered as hazardous substance.                                                                                                                               | 29 CFR 1910.1200<br>Occupational Safety and Health Standards                                                                          |
|   | This way up        | This is the correct upright position of the distribution packages for transport and/or storage.                                                                                          | ISO 780:2015(E)<br>Packaging-Distribution packaging-<br>Graphical symbols for handling and storage of packages                        |

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detax GmbH
Rapid Shape GmbH