**rs
VIVO**

RS VIVO IBT

rapidshape

bg	Инструкция за употреба	4
cs	Návod k použití	7
da	Brugsanvisning	10
de	Gebrauchsanweisung	13
el	Οδηγίες χρήσης	16
en	Instructions for use	19
es	Modo de empleo	22
et	Kasutusjuhend	25
fi	Käyttöohje	28
fr	Mode d'emploi	31
hr	Upute za uporabu	34
hu	Használati útmutató	37
it	Istruzioni per l'uso	40
lt	Naudojimo instrukcija	43
lv	Lietošanas instrukcija	46
nl	Gebruiksaanwijzing	49
no	Bruksanvisning	52
pl	Instrukcja stosowania	55
pt	Manual de instruções	58
ro	Instrucțiuni de utilizare	61
ru	Инструкция по применению	64
sk	Návod na použitie	67
sl	Navodila za uporabo	70
sv	Bruksanvisning	73
tr	Kullanım kılavuzu	76

ORDERING INFORMATION

RS VIVO IBT, 1000 g

Part no.

Standard

04722 [Rapid Shape part no.: RS006197]

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Resin for dental 3D printing

Indication

Orthodontic transfer trays for positioning brackets, bleaching trays

Patient target group

Persons being treated in the context of a dental procedure.

Intended users

Orthodontist, dentist, dental technician

Suitable for the following DLP printers/cleaning/post-exposure

see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Minimum material thickness for the design is 1.5 mm
- Processing temperature $23\text{ °C} \pm 2\text{ °C} / 73\text{ °F} \pm 4\text{ °F}$.

Safety information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Read and understand the safety data sheet!
- For users and/or patients

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- RS VIVO IBT is to be stored dry (at 15 °C - 28 °C / 59 °F - 82 °F) and protected from light. Minimal influence of light can already induce polymerisation.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxide. Some ingredients of RS VIVO IBT may cause allergic reactions in predisposed persons. In such cases refrain from using the product. RS VIVO IBT must only be inserted intraorally in completely polymerised state.

Adverse effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage

15 °C
59 °F

28 °C
82 °F

Processing

at 23 °C ± 2 °C / 73 °F ± 4 °F

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer.

① Construction process

Generation of a Print Job complying with machine and material parameters.

② Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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detax GmbH
Rapid Shape GmbH