

# D20+ cartridge

**The professional 3D printing solution for small labs with minimal investment.**

The D20+ cartridge features a large building area combined with the highest quality HD projection system. No compromise in quality or speed. With the cartridge, all consumables for 3D printing are provided, such as resin and the resin reservoir. Producing up to 7 parts per day, the D20+ cartridge is the ultimate deal for you to start professional 3D printing. Upgrading to a D20+ with an open material system is available at any time.

Performance parameters	D20+ cartridge
Building area	133 × 75 mm
Native pixel	+/-34 μm
Max. part height	90 mm
Light source	385 nm UV LED
Resolution	HD 1920 × 1080 px
Dimension (W × H × D)	335 × 541 × 349 mm
Connections	WLAN, TCP/IP, USB
Control	7" LCD-Display, touch control
Consumables	Cartridge, multi-use with 100 or 200 gr, resin and matching resin reservoir

Quick & easy multi-use cartridge system

Large printing area

Optional Force Feedback

Material identification (RFID)

Certified auto calibration (ACCS)

Remote control by TeamViewer®



Taylor made for small labs



Intelligent connectivity to RS wash, RS cure & cloud solutions



Clean & environmentally friendly disposal



Cartridges for many indications